


High Contrast Femtosecond Laser-driven K-shell X-ray Source in Imaging Application

Liming Chen (陈黎明), M. Kando, M. H. Xu, Y. T. Li, J. Koga,
M. Chen, H. Xu, X. H. Yuan, Q. L. Dong, Z. M. Sheng, S.V.
Bulanov, Y. Kato, T. Tajima, J. Zhang

- 1. Beijing National Laboratory of Condensed Matter Physics, Institute of Physics,
CAS, Beijing 100080, China*
- 2. Advanced Photon Research Center, Kansai Photon Science Institute, JAEA,
Kyoto 619-0215, Japan*
- 3. Department of Physics, Shanghai Jiao Tong University, Shanghai 200240,
China*



ICUIL 2008, Shanghai-Tongli




Outline

1. Current situation of fs laser-based hard x-ray source in medical imaging
2. Improved source using solid target
3. Improved source using cluster target
4. Summary

X-ray generation by ultrafast laser

- continuous bremsstrahlung emission (depends on temperature of hot electrons)
- characteristic x-ray emission lines (depend on target chemical composition)
- x-ray burst duration $\tau < 1$ ps (depends on optical pulse duration)
- x-ray focal spot size (depends on optical spot size and pulse duration); could be as small as ~ 2 μm

Source compare..



LB x-ray source


➤ The peak brightness is close to the 3rd generation synchrotron

With:

ps pulse duration

Monochromatic

Micro-m source size



Potential application

LB source is very useful for various fields:

Life Science: Atomic structure analysis of protein macromolecules; Mechanism of time-dependent biological reactions; Dynamics of muscle fibers

Materials Science: Precise electron distribution in novel inorganic crystals; Structural phase transition at high pressure / high or low temperature; Atomic and electronic structure of advanced materials of high T_c superconductors, highly correlated electron systems and magnetic substances; Local atomic structure of amorphous solids, liquids and melts

Chemical Science: Dynamic behaviors of catalytic reactions; X-ray photochemical process at surface; Atomic and molecular spectroscopy; Analysis of ultra-trace elements and their chemical states

Earth and Planetary Science: *In situ* X-ray observation of phase transformation of earth materials at high pressure and temperature; Structure of meteorites and interplanetary dusts

Environmental Science: Analysis of toxic heavy atoms contained in bio-materials; Development of novel catalysts for purifying pollutants in exhaust gases; Development of high quality batteries and hydrogen storage alloys



Industrial Application: Characterization of microelectronic devices and nanometer-scale quantum devices; Analysis of chemical composition and chemical state of trace elements; X-ray imaging of materials; Residual stress analysis of industrial products;

Medical Application: Application of high spatial resolution imaging techniques to medical diagnosis of cancers and capillaries

Difficulties for LB x-ray... (I) $K\alpha$ efficiency (η_K)

Using normal contrast laser...

Obstacles of $K\alpha$ x-ray enhancement when pulse width < 100 fs



η_K drops as higher laser intensities
---X-ray re-absorption

e.g: Ch. Reith et al., PRL 84, 4846(2000)


$L=0.3$, RA, Ag $I_{opt} \sim 10^{17}$ W/cm²

D. C. Eder et al., Appl. Phys. B 70, 211(1999)

η_K drops as shorter pulse width

e.g: M. Schnurer, et al., PRE, 61, 4394(2000);


J. Appl. Phys. 80, 5604(1996)



L. M. Chen, PoP 11, 4439(2004)

Saturation appears!!


Not optimized!!
10 times lower than theoretical prediction





D. Salzmann, PRE 65, 036402(2002)

Difficulties of LB x-ray... (II) Spectrum contrast


(c) spectrum after filter




(d) Image contrast on detector



(b) source spectrum



Example:
Experimental
mammography




Phantom 2 cm
25/75 (gland/adipose)
composition

• Imaging subject contrast becomes worse
When x-ray energetic tail increasing
---K α source with less background is necessary

	790 nm laser	X-ray tube
KT _h (keV)	25	28 kV _p
Imaging contrast	0.24 Worse !!	0.60


Spectrum contrast...

(a) Solid target




PoP 11, 4439(2004)

(b) Cluster target



PoP 11, 3491(2004)

(c) liquid micro-droplet





J. Phys. B 34, L313(2001)

	Solid target	droplet	clusters
Source C.	bad	worse	bad
K efficiency	$\sim 10^{-5}$ - 10^{-4}	$\sim 10^{-5}$	$\sim 10^{-6}$

LB x-ray source is seriously limited in imaging application !!
Stronger, cleaner and smaller source is necessary


How to improve?

Laser longitudinal compress via frequency doubling



X-ray source size detected vs. laser contrast
Smaller source size


Temporal profile of laser pulse




	800 nm	400 nm
Laser contrast ns	6×10^4	$> 10^{10}$
Laser contrast ps	1×10^4	$> 10^8$
Density gradient length (λ)	> 3	< 0.1
KT_h (keV)	> 150	~ 30
X-ray size (μm)	~ 100	~ 15

Effect of laser contrast ($I \sim 1 \times 10^{18} \text{ W/cm}^2$)

Overcomes the saturation limit...




Comparison for K α photon emission using 400 nm and 800 nm laser irradiation



K α efficiency depend as a function of laser intensity

- High contrast fs laser is more effective for K α photon generation
- K α efficiency does not show saturation at relativistic laser intensity using 400nm laser




Source contrast optimization...



X-ray spectrum obtained under
400 nm laser irradiation

- The energy cutoff show in x-ray spectrum of high contrast laser
----energetic x-ray tail suppressed

Vacuum Heating may be important mechanism





Calculation based on VH is more suitable for experiment than using RA

- Electron quiver amplitude: $X_{osc}/\lambda > L/\lambda$
- Energy cutoff (~ 20 keV) \Rightarrow prove VH
- quiver energy: $E_q = eE_L X_{osc} = 15$ keV
 \Rightarrow expel $J \times B$, RA (~ 100 keV)
- skin depth: $l_s \sim 10 v_{th}/v$, where $v_{th} = (2kT/m)^{1/2}$
 \Rightarrow expel SIB, ASE, STA ($l_s \ll v_{th}/v$)

\Rightarrow VH stimulated using high contrast laser


Control of hard x-ray flux via laser pulse shaping

Low contrast laser



KDP conversion efficiency vs. pulse duration under different pulse shape
 ---To get the right grating zero

High contrast laser




Kα conversion efficiency vs. pulse duration in differently skewed pulse shape


- η_K show significant chirp dependence
- The maximal $\eta_K \sim 4 \times 10^{-4}$, Similar to theoretical calculation limit.

Control of x-ray emission via laser pulse shaping

Calculation using HYADES



Temporal profile in different pulse duration/chirp



Plasma scale length in various pulse width and constant pulse energy
(laser intensity 1×10^{18} W/cm² on 60 fs)


$L(-100\text{fs}) = 0.1\lambda$ (optimal for VH)

$L(+100\text{fs}) = L(60\text{fs}) = 0.05\lambda$

The pre-plasma scale length can be controlled by the rising edge of high contrast laser pulse

Control of x-ray emission via pre-plasma

The peaked electron spectrum driven by high contrast laser




Integrated electron energy/absorption vs. density gradient


P-polarized laser with 30 optical cycles
at $1 \times 10^{18} \text{ W/cm}^2$ incident on target at 45°

- The level of VH stimulation is critically depend on L
- Preciously control of VH is possible via tuning the high contrast laser rising edge

Polarization dependence of hot electron generation




Angular distributions



- Electron show wider angular distribution along target normal. (evident of VH)
- Distinct target surface emission (P-pol.)
- Electron spectrum peaked.
- The second electron group appears (P-pol.) in spectrum, shows monoenergetic structure.

Spectra obtained in 125°



need more pure source
for imaging?





K-shell keV Ar X-ray source

Using TW laser

(Laser contrast 10^8)



High contrast Ar K-shell keV source





High contrast: peak/bkgd > 20
 $E_K/E_{total} > 70\%$
 Quasi-monochromatic
 Intense: $B = 10^{20}$ phs/s/mm²/mrad²
 Coherent (small size = 12 μm)
 Directional: $I_{Long.}/I_{Trans.} > 10$

■ Continuum background
 $\Sigma N_{K\alpha} / \Sigma N_{xray} \sim 10\%$
 -- polychromatic

Source optimization...




Long guiding (box 1.4x0.4 mm)




➤ Source flux enhanced when long channel formed

Phase-Contrast imaging




Sharp imaging...



SOD=38.5 cm
SDD=132 cm



Photon density on IP
 $\sim 5 \times 10^6 / \text{cm}^2$



Imaging edge greatly enhanced

New progress...

Imaging on CCD



Comparison ... **Ar keV flux x10 !!**

	Jlite-X	J-Karen
E_L (mJ)	200	800
τ (fs)	70	28
Contrast (to ps ASE)	10^6	$>10^9$ (adjustable)
I_{peak} (W/cm ²)	4×10^{17}	8×10^{18}
K α flux (phs/mrad ²)	1.2×10^3	1.1×10^4
	(phs/s/mm ² /mrad ²)	

Single-shot-imaging available !!

Photon flux on IP: 5×10^6 phs/cm²

Summary...

➤ The laser based hard x-ray source is optimized by using high contrast pulse:

Enhanced and controllable x-ray yield; Smaller x-ray source size

Suppressed energetic x-ray tail and energy cutoff

X-ray emission in shorter time scale (assumed)

Quasi-monoenergetic with low background

➤ In the field of *fs* laser physics, VH would be the main mechanism for hot electron generation. Precisely control the level of VH stimulation succeed via tuning the high contrast laser rising edge

➤ Extremely strong Ar K-shell x-ray source is archived. High contrast laser is the right choice to generate this source

Optimized x-ray source may allow imaging with...

- Improved specificity----higher contrast and higher spatial resolution
- Improved sensitivity----lower detection threshold
- Reduced dose on sample----better figure of merit (FOM)



Thank you